

OWEAR[®] OVERSHOES

Re-order codes & description:

OSN001	- Non-Skid	OSW001	- Waterproof CPE construction
OSN002	- Non Skid Large	OSW003	- Waterproof Over boot
OSR002	- Regular	OSW004	- Waterproof below Knee Over boot

Construction:

OSN001, OSN002 & OSR002

- Material*
- 40 gsm spunbond polypropylene with elastic opening
 - Anti-static and does not shed lint
 - Free of all foreign matter
 - OSN001 and OSN002 have an Isoprene non skid pattern on the sole
 - Latex free

OSW001

- Material*
- Polyethylene with elastic opening
 - Free of all foreign matter
 - Does not shed lint
 - Waterproof to fluid penetration

OSW003 & OSW004

- Material*
- 30gsm Polypropylene bonded to Polyethylene with elastic opening
 - Isoprene non-skid pattern on the sole
 - Antistatic
 - Free of all foreign matter
 - Waterproof to fluid penetration

OSW003 - Sonically bonded

OSW004 - Sonically bonded with stitching at ankle and calf

Specifications:

Colour – Blue

OSN001

Bottom length stretched	- 38cm min
Opening un-stretched	- 12cm min
Opening stretched	- 25cm min
Height	- 16cm min

OSN002

Bottom length stretched	- 45 min
Opening un-stretched	- 14cm min
Opening stretched	- 28cm min
Height	- 20.5cm min

OSR002

Bottom length stretched	- 38cm min
Opening un-stretched	- 12cm min
Opening stretched	- 25cm min
Height	- 16cm min

OSW001

Bottom length	- 35cm min
Opening un-stretched	- 15cm min
Opening stretched	- 35cm min
Height	- 14cm min

OSW003

Bottom length	- 39cm min
Opening un-stretched min	- 15cm min
Opening stretched	- 35cm min
Height	- 30cm min

OSW004

Bottom length	- 42cm min
Ankle opening un-stretched	- 15cm
Ankle opening stretched	- 24cm min
Calf opening un-stretched	- 11cm min
Calf opening stretched	- 24cm min
Height (sole to calf)	- 52cm min



Other - Latexfree
Inspection Tolerance - +5% variance
Shelf Life - 5 years

Compliance Standard:

- Class I
- ARTG Number - 245634

Recommended Storage:

Store in a clean, dry environment away from direct sunlight.

Packaging:

Primary (Selling Unit): 100 Owear[®] Overshoes in a dispenser

Shipper:

OSW001 – 10 selling units per shipper

Others – 6 selling units per shipper

Intended Purpose:

Owear[®] Overshoes are intended to protect and maintain hygiene and are therefore used in various industries such as hospitals, medical facilities and clean rooms. Being non-skid, they provide excellent grip to the feet while the waterproof overshoes provide fluid protection. Periodically check sole of overshoe to ensure abrasion has not reduced non-skid pattern (except OSR002) as this may reduce its effectiveness. Effectiveness of non-skid sole is also reduced in high wet areas eg showers, bathrooms. Caution and extreme care is advised in these areas.